


TREATREADY: ULTRA CONCENTRATED FEED SYSTEM



**STOP
BUYING
DRUMS**

PRODUCT APPLICATION GUIDE

- The **TreatREADY-1G** Ultra Concentrated System is a cost effective way to dispense ultra concentrates!
- **TreatREADY-1G** Ultra Concentrated System: Simply hang the TreatREADY Ultra Concentrated System on the wall, (hardware provided).
- Remove the child proof cap from the Ultra Concentrate
- Attach the blue S.A.F.E. cap to the top of the bottle, tighten snugly.
- Slip the bottle into the metal ring.
- Attach pump tubing to the top of the S.A.F.E. Cap and you're ready to go!
- The **TreatREADY-1G** Ultra Concentrated System, a self-contained, liquid Ultra Concentrate dispenser where there is no pouring or handing of the actual chemical treatment!
- The **TreatREADY-1G** Ultra Concentrated System can be used with Boiler Products: BoilerOxy-372, BoilerAid-383, BoilerAid-385, BoilerAid-386, BoilerAid-398, and BoilerAid-399.
- The **TreatREADY-1G** Ultra Concentrated System can be used with Cooling Products: KoolTower-121, KoolTower-123, KoolTower-125, KoolTower-130, and KoolLoop-210.
- The **TreatREADY-1G** Ultra Concentrated System can be used with microbiocides as well.

Application & Control

- Our Ultra Concentrated Products are packaged as 2 or 4 gallon cases, which are easily transported with grab & go handles and weighs lighter! This packaging is especially useful for the smaller account base where you may not wish to have excess inventory sitting around. If you have a large system, simply purchase 2 each of the 2 or 4 gallon cases! You can easily moderate inventory loads depending on application usage!